


User Manual 2
for
BL-18 B (Indian Beam line) at Photon Factory
Alignment of the 2nd Goniometer

Instructions for the users

Never allow the direct beam without proper absorber to fall on any type of detector.

Absorbers currently in use

- 64 - (64 x 20 micron Al foil)
- 32 - (32 x 20 micron Al foil)
- 16 - (16 x 20 micron Al foil)
- 8 - (8 x 20 micron Al foil)



Motors used in the liquid spectrometer:

All the goniometers at BL-18B are configured according to the right handed Cartesian co-ordinate system, where + y-axis represents the incoming beam direction.

*For changing incident energy, please refer **User Manual-1**.*


1st Goniometer motors







(Steering of Ge-crystal, will be mounted over this goniometer)

- filt - Filter before the ion chamber.
- th - Crystal rotation around the x-axis.
- tth - Detector rotation around the x-axis.
- alf - Crystal rotation around the beam direction, about y-axis, normally 'chi' in fourc.
- phi - Crystal in-plane rotation, rotation in the goniometer axis.
- mu - Whole goniometer rotation about the goniometer axis.
- s1slit - all the motors of slit-1 (primary beam slit, including the pseudo-motors for gap and center of the slits).
- s2slit - all the motors of slit-2 (diffracted beam slit, including the pseudo-motors for gap and center of the slits).
- tbx - Goniometer x-motion motor.
- tbz - Goniometer z-motion motor.
- samx - Sample x-motion motor.
- samy - Sample y-motion motor.
- samz - Sample z-motion motor.
- bm_h - Beam monitor horizontal.
- bm_v - Beam monitor vertical.

2nd Goniometer motors (Sample will be mounted over this goniometer)

sh –	Whole goniometer height motion motor.
sx –	Sample x-motion motor.
sy –	Sample y-motion motor.
sz –	Sample z-motion motor (for fine tuning, this motor is sensitive to the load on the sample stage).
sth –	Sample in-plane motion around the z-axis.
beta –	Detector rotation around x-axis.
stth –	Detector in-plane rotation around the z-axis.
xtilt –	Sample tilt around x-axis.
ytilt –	Sample tilt around y-axis.
ih –	height of flight path.
ir –	flight path rotation around x-axis.
ome –	Vertical sample stage horizontal motion.
hgonio –	Air pad motion for total goniometer motion.
is –	slit assembly between the diffractometers (horizontal slit was connected).





```

230. SURF> wa
Current Positions (user, dial)
SI_left Inc_beam Sample_he SI_left 2Theta Det_hor_r Liq_hor_r Phi
sll ih sh sll tth stth sth phi
0.4086 -189.8616 0.0000 0.4086 0.0000 0.5000 0.0000 0.0000
-0.5964 -3.5786 130.2227 -0.5964 0.0000 -0.1380 0.0000 -5.7366
Inc_beam SI_left SI_right SI_bottom SI_top S2_bottom S2_left S2_right
ir sll slr slb slt s2b s2l s2r
-0.0195 0.4086 1.5914 -1.2386 1.3886 0.1000 1.7019 0.2981
0.0077 -0.5964 0.5664 0.2819 -0.3269 0.1608 1.0519 -0.1118
S2_top Slv_gap Slv_cen Slh_gap Slh_cen S2v_gap S2v_cen S2h_gap
s2t slvgap slvcen slhgap slhcen s2vgap s2vcen s2hgap
0.1000 0.1500 1.3136 2.0000 -0.5914 0.2000 0.0000 2.0000
6.2219 0.1500 1.3136 2.0000 -0.5914 0.2000 0.0000 2.0000
S2h_cen Table_X Table_Z Mu Alpha Theta Sample_X Sample_Y
s2hcen tb_x tb_z alf th samx samy
0.7019 -0.2186 -1.4500 0.0000 -90.0000 0.2000 18.0000 0.0000
0.7019 500.0000 -0.2500 0.0000 -0.0323 0.1341 18.0000 0.0000
Sample_Z Beam_mon Beam_mon Filter Liq_sam_x Liq_sam_y Liq_sam_z Liq_ytilt
samz bm_h bm_v filt sx sy sz ytilt
17.3520 0.0000 0.0000 0.0000 0.0000 4.0000 13.0150 0.0000
11.9620 0.0000 0.0000 0.0000 0.0000 4.0000 13.0150 0.0000
Liq_xtilt Det_ver_r Air_wheel Ver_stage Inc_slit Inc_slit
xtilt beta hgonio ome isb ist
-1.0000 0.0000 -0.1950 -8.2660 0.5000 0.5000
-1.0000 -13.7543 -0.1950 8.2660 -12.6260 -4.3000

```

Use “*surf*” for this alignment.

To align the 1st goniometer, w.r.t direct beam, please refer user manual-1.

Mount the Ge(111) crystal and find the diffraction beam (keeping absorbers ~ 12° 2-theta for 18keV energy) and do the fine adjustment using ‘*th*’, ‘*tth*’, ‘*alf*’, and ‘*phi*’ (‘*lup*’ command).

Move ‘*alf*’ (‘*chi*’) 90° and shift the detector from 1st goniometer to second goniometer.

Keep the 2theta arm at 90° so that it will not come in the path of 2nd goniometer.

Focus the detector to the center of the 2nd goniometer.

‘*samz*’ of the 1st goniometer should be moved down so that the sample is out of the direct beam in order to align the 2nd goniometer.

There are two pins on the 2nd goniometer one is on the vertical sample mounting stage and other is on the normal (horizontal) sample mount. Use “*ome*” command to bring vertical stage of the 2nd goniometer to the goniometer center and adjust height of the pin on the horizontal sample stage such that it comes in the direct beam then align the goniometer to the direct beam.

Adjust the height of the 2nd goniometer such that the direct beam passes through the detector and cutting the pin at the center of the goniometer (use all the absorbers).

Move ‘*samz*’ to the original position (half cut) and ‘*mu*’ to » 6° (in case of E = 18 keV) for diffraction peak. Use florescent screen to view the diffracted beam and fine tune ‘*mu*’.

Slowly move the 2nd goniometer to the diffraction peak using the air pad motor (‘*hgonio*’) and the florescent screen.

After finding the diffracted beam fine tune ‘*mu*’, ‘*phi*’ and ‘*alf*’ to get maximum intensity. Then align the 2nd goniometer to the diffracted beam.

For aligning the 2nd goniometer to the diffracted beam, fine tune ‘*alf*’ and ‘*beta*’ while reducing the detector slits and align ‘*stth*’ as well to the direct beam.

Keep the flight path, ion chamber and the slit assembly (manual and motor controlled) between the two goniometers and adjust to the diffracted beam.